

# ***BC-401F FitScan Travel Scale***

## ***Body Composition Monitor***



***Easily monitor your progress  
anywhere you go***

***Hard case is light weight but  
rugged for durability***

***10 full-body compositions  
measurements***

***Built in Bluetooth wireless  
allows easy connection to  
your smart devices***

***Includes Tanita's **FREE** Healthy  
Edge mobile app***

***Large easy to read display even  
in low light***



Contract Holder  
Contract 36F79719D0239



## ***BC-401F FitScan Travel Scale*** ***FULL Body Composition Monitor***

The BC-401F FitScan Travel Scale is the perfect tool for monitoring the effectiveness of any fitness or health program anywhere you go. This Full Body Composition Monitor includes a **hard travel case with purchase** and displays 10 body composition measurements. The extra-large LCD color display is easy-to-read read even in low light. The BC-401F FitScan Travel Scale can easily transfer your full body composition measurements to your personal smart device using Tanita's **FREE Healthy Edge mobile app**, since it is equipped with Bluetooth Wireless. Up to 4 users can record and review their results (each user must have own smart device). This travel scale is **FDA cleared** and **has our medically trusted Bioelectrical Impedance Analysis technology (BIA)**.



### **Measurements Include:**

- ✓ Weight
- ✓ Body fat %
- ✓ Body water %
- ✓ Muscle mass
- ✓ Body mass index
- ✓ Metabolic age (12 - 90)
- ✓ Bone mass
- ✓ Visceral fat rating
- ✓ Basal metabolic rate
- ✓ Physique rating

### **Specifications:**

- ✓ Capacity: 330lb (150 kg)
- ✓ Weight increments: 0.2lb (0.1kg)
- ✓ Platform Size: 8.5" x 12.5"
- ✓ Bluetooth Wireless
- ✓ Easy-to-read LCD display
- ✓ 4 person memory
- ✓ Guest Mode & Athlete Mode
- ✓ Recall function
- ✓ Power Supply: 4 x AA batteries (*included*)
- ✓ 3 year warranty