

BC-533

BODY COMPOSITION MONITOR



Popular Tanita InnerScan monitor

9 full-body measurements

DCI calculation for weight loss

Rugged tempered glass platform

Distinctive round InnerScan design

4 person memory

The Tanita BC-533 InnerScan Body Composition Monitor uses FDA cleared Bioelectric Impedance Analysis (BIA) technology in an innovative, circular design, which was designed for people with distinct taste and a passion for personal health. It delivers hi-tech readings in categories like weight, body fat percentage, body water percentage, muscle mass, bone mass, and more. BIA is quick and non-invasive, and is one of the most thorough and reliable ways to measure body composition, clinically comparable to DXA (dual-energy x-ray absorptiometry) and hydrostatic (underwater) weighing.

The BC-533 is one of Tanita's bestselling body composition monitors, and is the perfect companion to help you achieve your health and wellness goals. The platform has a weight capacity of 330 pounds, and operates for greater than 2 years on four AA batteries (included).

Measurements Include:

- ✓ Weight
- ✓ Body water %
- ✓ Body fat %
- ✓ Muscle mass
- ✓ Physique rating
- ✓ Metabolic age
- ✓ Bone mass
- ✓ Visceral fat rating
- ✓ Daily calorie intake (DCI)

Specifications:

- ✓ Capacity: 330lb (150 kg)
- ✓ Weight increments: 0.2lb (0.1kg)
- ✓ Body Fat Increments: 0.1%
- ✓ Tempered glass platform
- ✓ Metal electrodes
- ✓ Recall function
- ✓ Memory for up to 4 users
- ✓ Weight only function
- ✓ Guest Mode & Athlete mode
- ✓ Power Supply: 4 x AA (included)
- ✓ Dimensions: 13.4" x 13.4" x 1.3"
- ✓ Weight: 5.7 lb