

# BC-568

## SEGMENTAL BODY COMPOSITION MONITOR



**Full Body Composition (weight, body fat%, body water%, muscle mass, physique rating, and more...)**

**Segmental Body Composition (Right Arm, Left Arm, Right Leg, Left Leg, and trunk)**

**Automatic scrolling of results: weight, body fat%, and segmental body fat%**

**Recall function: weight, body fat%, body water%, muscle mass, BMR, bone mass, and visceral fat**

**Graph function: body fat%, body water%, muscle mass, BMR, bone mass, and visceral fat**

The BC-568 Segmental Body Composition Analyzer from Tanita represents the most advanced home health monitoring product, a whole new level in measuring your Competitive Edge. Technology found in the BC-568 is based on the same highly accurate, FDA cleared, advanced technology found in Tanita's professional line of segmental body composition analyzers. Utilizing Tanita's exclusive bioelectric impedance analysis (BIA) technology this monitor not only measures weight and body fat, it also tells you your muscle mass, body water %, basal metabolic rate (BMR), metabolic age, bone mass, visceral fat and rates your physique!

The BC-568 uses innovative retractable handgrip electrodes, along with standard feet electrodes. These eight BIA electrodes allow the BC-568 to go well beyond standard body composition; giving you a segmental review of what you're made of. Why segmental? Simply put, segmental gives you more information than standard body composition monitors. This totally unique product gives individual body composition readings for each body segment- trunk, right arm, left arm, right leg and left leg. The BC-568 is especially useful for anyone who is monitoring the balance of left and right side of body or trying to build or rehabilitate a particular part of one's body.

### Measurements Include:

- ✓ Weight
- ✓ Body water %
- ✓ Body fat % **(total & segmental)**
- ✓ Muscle mass **(total & segmental)**
- ✓ Physique rating
- ✓ Metabolic age
- ✓ Bone mass
- ✓ Visceral fat rating
- ✓ Daily calorie intake (DCI)

### Specifications:

- ✓ Capacity: 330lb (150 kg)
- ✓ Weight increments: 0.2lb (0.1kg)
- ✓ Body Fat Increments: 0.1%
- ✓ Recall function: weight, body fat %, body water %, muscle mass, DCI, bone mass and visceral fat
- ✓ Graph function for: body fat % (total and segmental), body water %, bone mass, visceral fat and DCI
- ✓ Graph recall: calendar days / 52 weeks / 36 months
- ✓ Guest Mode
- ✓ Backlit icon buttons
- ✓ Power Supply: 4 x AA *(included)*
- ✓ Dimensions: 13.4" x 14.4" x 1.7"
- ✓ Weight: 6.6 lb