



BODY COMPOSITION READINGS

- Weight
- Body fat %
- Body water %
- Muscle mass
- Physique rating
- Daily Caloric Intake
- Metabolic age
- Bone mass
- Visceral fat

BC-585F BODY COMPOSITION MONITOR

Bring out your inner athlete with the BC-585F Body Composition Monitor. The BC-585F features an innovative design while giving you access to all the best information through BIA Technology: weight, body fat percentage, body water percentage, muscle mass, bone mass, and so much more. This is a great fitness tool to keep on display in your bathroom - it's as fashionable as it is functional - and one you can proudly display because of its aesthetic appeal.

Specifications:

- Capacity: 330 lb (150 kg)
- Weight Increments: 0.2 lb (0.1 kg)
- Body Fat Increments: 0.1 %
- Easy-to-read 2.25" LCD display
- Tempered glass platform
- 4 person memory
- Guest Mode & Athlete Mode
- Recall function: weight, body fat %, body water %, and muscle mass
- Power: 4 x AA batteries (included)
- 1 year warranty
- Dimensions: 12.6" x 11.8" x 1.7"
- Weight: 5.9 lb