

# **BC-601FS**

## **SEGMENTAL BODY COMPOSITION MONITOR**



**Entry-level segmental body composition analyzer**

**Measures arms/legs/trunk fat and lean tissue**

**19 measurements total (full-body and segmental)**

**Monitor body changes and balance**

**Innovative display in handle for easy viewing of results**

***SD Card  
Capability***

## **BC-601FS**

### **SEGMENTAL BODY COMPOSITION MONITOR**

The FitScan BC-601FS is the perfect tool for monitoring the effectiveness of any fitness or health program. This 8 electrode Body Composition Monitor displays 19 readings; 9 whole body AND 10 segmental (arms, legs and trunk area) so you can see exactly how your body is responding to specific training and lifestyle changes. From tracking body fat and muscle mass trends, to recommending calorie intake and indicating hydration levels this Monitor helps you tailor your exercise and dietary requirements to achieve your goals. The extra large LCD color display is easy-to-read and the retractable cord ensures the handset fits snugly into the main unit for easy storage. Includes SD Card to record all measurements for further analysis. Up to 4 users can record and review their results, which are stored for 6 months.



#### **Measurements Include:**

- ✓ Weight
- ✓ Body fat %
- ✓ Body water %
- ✓ Muscle mass
- ✓ Daily Caloric Intake
- ✓ Metabolic age
- ✓ Bone mass
- ✓ Visceral fat rating
- ✓ Segmental body fat %
- ✓ Segmental muscle mass
- ✓ Body mass index

#### **Specifications:**

- ✓ Capacity: 330lb (150 kg)
- ✓ Weight increments: 0.2lb (0.1kg)
- ✓ Body Fat Increments: 0.1%
- ✓ Easy-to-read LCD display
- ✓ 4 person memory
- ✓ Guest Mode & Athlete Mode
- ✓ Recall function
- ✓ Power Supply: 4 x AA batteries (*included*)
- ✓ 90 day warranty