

# BC-730F

## BODY COMPOSITION MONITOR



**Compact lightweight design**

**9 full-body measurements**

**Auto user recognition**

**5 person memory**

**Comes in 4 colors**



The BC-730F is super compact and lightweight with an easy-to-read display. It has an auto-recognition feature so after you have programmed your details, the monitor will recognize you each time you step on, shortening the measurement time to a few seconds. Providing 9 full-body measurements, the BC-730F is the perfect monitor for those with limited space, as well as a travel companion for those who want to track their body composition away from home.

### Measurements Include:

- ✓ Weight
- ✓ Body fat %
- ✓ Body water %
- ✓ Muscle mass
- ✓ Basal metabolic rate (BMR)
- ✓ Metabolic age
- ✓ Bone mass
- ✓ Visceral fat rating
- ✓ Body Mass index (BMI)

### Specifications:

- ✓ Capacity: 330lb (150 kg)
- ✓ Weight increments: 0.2lb (0.1kg)
- ✓ Body Fat Increments: 0.1%
- ✓ Easy-to-read 2 line
- ✓ 5 person memory with recall
- ✓ Guest Mode & Athlete Mode
- ✓ Power Supply: 4 x AA batteries (*included*)
- ✓ 1 year warranty
- ✓ Dimensions: 8.5" x 10.2" x 1.3"
- ✓ Weight: 2.0 lb