

BF-684W

Body Fat Monitor



Tanita's best body fat / body water scale

Highest capacity of 440 pounds

Oversized tempered glass platform

Beautiful blue backlit LCD and icons

Includes special mode for athletes

Body fat monitoring is an essential tool for precisely measuring your body's inner health. Tanita Body Fat Monitors provide vital in-depth information, that allows you to accurately observe the impact of a diet and fitness program.

Tanita's new development allows users to monitor body composition on clear surfaces without visible electrodes. This innovation will add to the aesthetic beauty of Tanita's unique range of tempered glass designs while providing the same high level of accuracy. This transparent electrode scale uses Bioelectrical Impedance Analysis (BIA) technology to monitor your body fat and body water. Model features higher weight capacity 440lbs, wider platform size 13.4" and blue backlit LCD display.

Measurements:

- ✓ Weight
- ✓ Body water %
- ✓ Body fat %

Specifications:

- ✓ Capacity: 440lb (200 kg)
- ✓ Weight increments: 0.2lb (0.1kg)
- ✓ Body Fat Increments: 0.1%
- ✓ Tempered glass platform
- ✓ Transparent electrodes
- ✓ Blue backlight display
- ✓ Memory for up to 2 users
- ✓ Guest Mode & Athlete Mode
- ✓ Recall function
- ✓ LED icon guide
- ✓ Power Supply: 4 x AA (included)
- ✓ Dimensions: 12.6" x 13.4" x 1.7"
- ✓ Weight: 5.9 lb

New Blue Backlight Display Feature



Size: 2.5" w x 2.3" h