

# **IRONMAN® DC-240IM**

## *PROFESSIONAL BODY COMPOSITION ANALYZER*



*The first professional grade  
**IRONMAN** body composition  
monitor, **THE STANDARD**  
by which all others are measured*

*Bundled with Tanita's award  
winning Health Ware® Software,  
**MEASURE - GRAPH - ANALYZE***

*The perfect tool to monitor health  
and wellness as well as athletes  
training for an event, or recovering  
from injury*

*Monitor healthy weight by analyzing  
fat loss versus muscle/water loss*

*The most cost effective system  
providing professional results,  
from the most trusted brand, Tanita*

### DC-240IM

#### PROFESSIONAL BODY COMPOSITION ANALYZER

The **IRONMAN** brand is the standard by which all other products are measured, only the premier products are allowed to carry the **IRONMAN** logo. The DC-240IM is a bundled package that provides everything needed to professionally monitor health and fitness. Included is Tanita's award winning Health Ware® software, and a USB cable to connect to your Windows based PC (single license).

The DC-240IM is an FDA cleared, lightweight, highly portable, body composition monitor. Weighing just 10 pounds with a built-in handle, the platform easily accommodates individuals up to 440 pounds. Tanita Health Ware® software provides custom data recording, tracking and analysis for an unlimited number of users. Reports are easily created using standard color printers, no special forms are required. The DC-240IM and Tanita Health Ware® combination provides detailed analysis and long term trend monitoring of body composition data, vital in assessing the effectiveness of any health or fitness program.

#### Features:

Data Shown on the LCD Display:

- ✓ Weight
- ✓ Body Fat %
- ✓ Body Water %
- ✓ BMI


Additional data available using Tanita Health Ware®:

- ✓ Fat mass
- ✓ Fat free mass
- ✓ Body water mass
- ✓ Muscle mass
- ✓ Bone mass
- ✓ Visceral fat rating
- ✓ Metabolic age
- ✓ Muscle score
- ✓ Daily caloric intake (DCI)
- ✓ Basal metabolic rate (BMR)
- ✓ Ideal body weight
- ✓ Degree of obesity
- ✓ BMR score
- ✓ Impedance

#### Specifications:

- ✓ Capacity: 440lb (200 kg)
- ✓ Weight Graduation: 0.2lb (100 g)
- ✓ Body Fat Increments: 0.1%
- ✓ USB Type B Output
- ✓ Dimensions: 17.2" x 13.4" x 2.1"
- ✓ Weight: 10.0lbs
- ✓ Power Supply:
  - AC Adapter (included)
  - or 6 x AA (not included)

#### Tanita Health Ware® detailed Body Composition Report



Optional Accessories:



C-300 Carrying Case