


WB-800Hplus

DIGITAL SCALE WITH BMI FUNCTION



***Hospitals, Clinics, Spas,
Gyms, Clubs, these
scales are the perfect
weighing solution***

***660 lb / 300kg weight
capacity suitable for all
patients***

***EMR ready with
straightforward digital
outputs***

***Step-on, activation, no
buttons to push***

***Low level stable
platform for elderly or
obese patients***

***Swivel easy-to-read
LCD displays allowing
confidential readings***

Integrated Stadiometer

WB-800Hplus

DIGITAL SCALE WITH BMI FUNCTION

The WB-800Hplus scale is perfect for those instances where absolute measurement of weight gain or loss is a necessity. This incredibly durable, yet portable, professional grade scale is perfect for any busy practice, gym or clinic. Features an extended capacity and **oversized low-profile platform**, with instant-on activation of digital weight and BMI to accommodate a wide range of clients and patients with fast, accurate results on a jumbo LED display. With an integrated stadiometer on column mounted display provides efficient measurement of height, suitable for both children and adults. Includes rotating display and platform castors. The RS-232 and USB outputs provide straight-forward digital data for EMR integration.

Specifications:

- ✓ **Capacity:** 660 lb (300 kg)
- ✓ **Weight Increments:** 0.2 lb (0.1 kg)
- ✓ **Height Range:**
25"–84" x 1/4" / 64–214 cm x 0.5 cm
- ✓ **Dimensions and Weight:**
WB-800Hplus: 23.4" x 15.6" x 56.1", Wt 25.0lb
- ✓ **Power:** AC adapter (included)
4 x AA batteries (not included)
- ✓ **Platform:** 15.3" x 15.6" x 1.9"
- ✓ **Display:** 2.4" x 3.9"
- ✓ **Output:** RS-232 and USB (type B)
- ✓ **Comes with a 3 year warranty**

Compatible With:

- ✓ Tanita Health Ware™ Software
- ✓ Tanita's WA-232 (accessory),
RS-232 to ANT+
Radio Wireless Adapter



Wirelessly transmit measurements to Tanita Health Ware™ software, then export to Excel or .csv data.

- ✓ RS232 Bluetooth Radio Wireless Adapter
Use for secure extended range up to 300 feet!



Plugs into the USB in your computer



Plugs into the RS-232 port on the scale

WB-800Hplus is also Available in kg-Only Mode