

WB-800Pplus

DIGITAL SCALE WITH BMI FUNCTION



***Hospitals, Clinics, Spas,
Gyms, Clubs, these
scales are the perfect
weighing solution***

***660 lb / 300kg weight
capacity suitable for all
patients***

***EMR ready with
straightforward digital
outputs***

***Step-on, activation, no
buttons to push***

***Low level stable
platform for elderly or
obese patients***

***Swivel easy-to-read
LCD displays allowing
confidential readings***

WB-800Pplus

DIGITAL SCALE WITH BMI FUNCTION

The WB-800Pplus scale is perfect for those instances where absolute measurement of weight gain or loss is a necessity. This incredibly durable, yet portable, professional grade scale is perfect for any busy practice, gym or clinic. Features an extended capacity and **oversized low-profile platform**, with instant-on activation of digital weight and BMI to accommodate a wide range of clients and patients with fast, accurate results on a jumbo LED display. Column mounted display rotates 90 degrees for ease of viewing in any situation. Castors on the base of the platform make moving from room to room effortless. The RS-232 and USB outputs provide straight-forward digital data for EMR integration.

Specifications:

- ✓ **Capacity:** 660 lb (300 kg)
- ✓ **Weight Increments:** 0.2 lb (0.1 kg)
- ✓ **Dimensions and Weight:**
WB-800Pplus: 23.8" x 15.6" x 47.1", Wt 22.2lb
- ✓ **Power:** AC adapter (included)
4 x AA batteries (not included)
- ✓ **Platform:** 15.3" x 15.6" x 1.9"
- ✓ **Display:** 2.4" x 3.9"
- ✓ **Output:** RS-232 and USB (type B)
- ✓ **Comes with a 3 year warranty**

Compatible With:

- ✓ Tanita Health Ware™ Software
- ✓ Tanita's WA-232 (accessory), RS-232 to ANT+ Radio Wireless Adapter



Wirelessly transmit measurements to Tanita Health Ware™ software, then export to Excel or .csv data.

- ✓ RS232 Bluetooth Radio Wireless Adapter
Use for secure extended range up to 300 feet!



Plugs into the USB in your computer



Plugs into the RS-232 port on the scale

WB-800Pplus is also Available in kg-Only Mode